

Test	Test Description	REFLEX TO	Test Description	NOTES
CBCWD	CBC with Automated Differential	MDIFF	Manual Performed Differential	Certain abnormal WBC, RBC, or platelete values
SYIGG	Treponema Pallidum Antibodies (IgG, IgM)	RPR	Rapid Plasma Reagin	If SYIGG and RPR results are discordant, RPR reflex to TPPA (reference send out; manually ordered and manually received)
HBSAG	Hepatitis Bs Antigen	HBSGC	Hepatitis Bs Ag, Confirmation	
HIV12	HIV 1/2 Ag/Ab	HIVAC	HIV 1,2 Antibody Confirmation	Patients less than 2 years of age or if the confirmation is discordant, reflex to HIV12 RNA Detect, LAB10757 (automatically ordered, manually received)
TSH	TSH	LAB127	FT4	Any L or H flag will reflex to FT4 (automatically ordered and received)
ANA	Anti-nuclear Antibody	ANA Titer, DSDNA, SSA, SSB, SM, SM/RNP	ANA Panel	
HCVAB	Hepatitis C Antibody	HCVQN	Hepatitis C RNA Detect/Quant	Automatically ordered, manually received, cannot be run on patients <18 years old
GHGB	Hemoglobin A1c	HBA1C	Hemoglobin A1c	Variant result, reflex to reference lab hemoglobin A1c
CDIFF	Clostridioides Difficile	CDPCR	C DIFF TOXIN BY PCR	Reflexes when CDIFF EIA kit is indeterminate
UMAC	Urine Macroscopic	UMIC	Urine Microscopic	Unit = Positive UPRO = 1+, 2+, or 3+ ULEU = Trace or 1+, 2+, or 3+ UBLD = 1+, 2+, or 3+ UAPP = Cloudy UGLU = 3+ or 4+
UAR	Urinalysis with reflex to culture	CXURN	Urine Culture	UWBC is >10 or >=100
HBABL	Hepatitis B Surface Antibody			Cannot be run on patients <30 days (neonates)
HAVG	Hepatitis A Antibody, IgG			Cannot be run on patients <4 years old; no reference reflex available